


**BIONET** core scaffolds -  
A Key step towards discovery



6ece\_Ter 2009  
Core Scaffolds Update

 **Key Organics**  
*The right side of the Equation*

**LD-0080**

NEW

GBP (£)


 $C_{13}H_{12}BrNO_3$ 

MW 310.1

Mp/Bp

224 - 225

Purity &gt;95%



500mg 136

1g 182

5g

10g

5-(benzyloxy)-3-bromo-2-(hydroxymethyl)-4(1H)-pyridinone

**LD-0048**

NEW

GBP (£)


 $C_{13}H_{11}BrINO_2$ 

MW 420.0

Mp/Bp

133 - 136

Purity &gt;95%



500mg 168

1g 229

5g

10g

[5-(benzyloxy)-6-bromo-4-iodo-2-pyridinyl]methanol

**MD-0033**

NEW

GBP (£)


 $C_{13}H_{12}INO_2$ 

MW 341.2

Mp/Bp

215 - 217

Purity &gt;95%



500mg 101

1g 137

5g

10g

3-(benzyloxy)-5-iodo-2-methyl-4(1H)-pyridinone

**MD-0047**

NEW

GBP (£)


 $C_9H_5BrClNO$ 

MW 258.5

Mp/Bp

260 - 261

Purity &gt;95%



500mg 202

1g

5g

10g

5-bromo-2-chloro-1H-indole-3-carbaldehyde

**FD-0056**

NEW

GBP (£)


 $C_9H_4BrClN_2O_2$ 

MW 287.5

Mp/Bp

182 - 183

Purity &gt;95%



500mg 101

1g 137

5g

10g

6-bromo-4-chloro-8-nitroquinoline

**LD-0097**

NEW

GBP (£)


 $C_{13}H_{14}BrNO$ 

MW 280.2

Mp/Bp

197 - 198

Purity &gt;95%



500mg 101

1g

5g

10g

5'-bromo-1',2'-dihydrospiro[cyclohexane-1,3'-indole]-2'-one

**MD-0046**

NEW


 $C_{11}H_{11}BrN_2O$ 

MW 267.1

Mp/Bp

205 - 206

Purity &gt;95%



GBP (£)

500mg 136

1g

5g

10g

6-bromo-2-(dimethylamino)-1H-indole-3-carbaldehyde

**MD-0055**

NEW


 $C_8H_6BrIO_3$ 

MW 356.9

Mp/Bp

197 - 199

Purity &gt;95%



GBP (£)

500mg 136

1g 182

5g 456

10g

5-bromo-3-iodo-2-methoxybenzenecarboxylic acid

**LD-0033**

NEW


 $C_7H_8BrNO_2$ 

MW 218.1

Mp/Bp

75 - 78

Purity &gt;95%



GBP (£)

500mg 101

1g 137

5g

10g

(6-bromo-5-methoxy-2-pyridinyl)methanol

**MD-0041**

NEW


 $C_9H_6BrF_3N_2O_3$ 

MW 327.1

Mp/Bp

101 - 103

Purity &gt;95%



GBP (£)

500mg 68

1g 92

5g

10g

N-[4-bromo-2-nitro-5-(trifluoromethyl)phenyl]acetamide

**MD-0250**

NEW


 $C_9H_{12}ClN_3O$ 

MW 213.7

Mp/Bp

152 - 153

Purity &gt;95%



GBP (£)

500mg 101

1g 137

5g

10g

2-chloro-N'-hydroxy-6-isopropyl-3-pyridinecarboximidamide

**MD-0239**

NEW


 $C_{12}H_{10}ClN_3O$ 

MW 247.7

Mp/Bp

182 - 183

Purity &gt;95%



GBP (£)

500mg 136

1g 182

5g

10g


2-chloro-N'-hydroxy-6-phenyl-3-pyridinecarboximidamide

**MD-0726**

NEW

GBP (£)

$C_9H_9ClN_2$   
 MW 180.6  
 Mp/Bp  
 25 - 26  
 Purity >95%



500mg 101  
 1g 137  
 5g 342  
 10g


2-chloro-6-isopropylnicotinonitrile

**MD-0715**

NEW

GBP (£)

$C_8H_5ClN_2O_4$   
 MW 228.6  
 Mp/Bp  
 283 - 285  
 Purity >95%  
 [116862-22-7]



500mg 50  
 1g 69  
 5g 173  
 10g


6-chloro-7-nitro-2H-1,4-benzoxazin-3(4H)-one

**MD-0729**

NEW

GBP (£)

$C_8H_7ClN_2O_3$   
 MW 214.6  
 Mp/Bp  
 204 - 206  
 Purity >95%



500mg 68  
 1g 92  
 5g  
 10g


6-chloro-7-nitro-3,4-dihydro-2H-1,4-benzoxazine

**LD-0734**

NEW

GBP (£)

$C_8H_8ClN_3O_4$   
 MW 245.6  
 Mp/Bp  
 122 - 125  
 Purity >95%



500mg 168  
 1g  
 5g  
 10g


ethyl 4-amino-6-chloro-5-nitronicotinate

**LD-0719**

NEW

GBP (£)

$C_7H_6FNO_3$   
 MW 171.1  
 Mp/Bp  
 69 - 70  
 Purity >95%  
 [63878-73-9]



500mg 33  
 1g 45  
 5g 113  
 10g 170


(2-fluoro-5-nitrophenyl)methanol

**MD-0704**

NEW

GBP (£)

$C_8H_{10}N_2O_3$   
 MW 182.2  
 Mp/Bp  
 157 - 159  
 Purity >95%



500mg 68  
 1g 92  
 5g  
 10g

[2-(methylamino)-5-nitrophenyl]methanol

**FD-0072**

NEW


 $C_{11}H_7BrClNO_2$ 

MW 300.5

Mp/Bp

226 - 228

Purity &gt;95%



GBP (£)

500mg 101

1g 137

5g

10g

methyl 6-bromo-4-chloro-8-quinolinecarboxylate

**LD-0027**

NEW

 $C_8H_6BrIO_3$ 


MW 356.9

Mp/Bp

126 - 128

Purity &gt;95%

[18071-51-7]



GBP (£)

500mg 68

1g 92

5g 229

10g

methyl 5-bromo-2-hydroxy-3-iodobenzenecarboxylate

**MD-0013**

NEW

 $C_8H_6BrNO_5$ 


MW 276.0

Mp/Bp

147 - 149

Purity &gt;95%

[91983-31-2]



GBP (£)

500mg 68

1g 92

5g 229

10g

methyl 5-bromo-2-hydroxy-3-nitrobenzenecarboxylate

**MD-0023**

NEW


 $C_9H_8BrIO_3$ 

MW 371.0

Mp/Bp

65 - 67

Purity &gt;95%



GBP (£)

500mg 101

1g 137

5g 342

10g

methyl 5-bromo-3-iodo-2-methoxybenzenecarboxylate

**MD-0238**

NEW


 $C_{16}H_{16}N_4$ 

MW 264.3

Mp/Bp

78 - 79

Purity &gt;95%



GBP (£)

500mg 136

1g 182

5g

10g

6-phenyl-2-piperazinonicotinonitrile

**MD-0732**

NEW


 $C_6H_{11}Cl_2N_3$ 

MW 196.1

Mp/Bp

260 - 264

Purity &gt;95%



GBP (£)

500mg 168

1g 229

5g 572

10g 857

2H,4H,5H,6H,7H-pyrazolo[4,3-c]pyridine dihydrochloride

**MD-0212**

**NEW**


C<sub>21</sub>H<sub>24</sub>N<sub>4</sub>O<sub>2</sub>

MW 364.5

Mp/Bp

158 - 159

Purity >95%



**GBP (£)**

**500mg 136**

**1g 182**

**5g**

**10g**

**tert-butyl 4-(3-cyano-6-phenyl-2-pyridinyl)tetrahydro-1(2H)-pyrazinecarboxylate**

---